



1000 Corporate Drive, Marshfield, WI 54449
1-844-484-4998 · www.argondentalusa.com


Argon K3Pro™ 4.0 - 4.5mm Rapid Precision Surgical Kit

Drill Usage Chart

IMPORTANT NOTE:

To maintain optimal performance, Argon Dental USA recommends replacing drills after 20 uses. Please use this chart as a tracking device. Check the box which corresponds to the drill, each time it is used.

4.0-4.5mm Punch, Starter and Trepan



	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	Product #	
4.0-4.5 Punch																											GS_RS4.0-4.5
4.0-4.5 Starter																											RS_PF4.0-4.5
4.0-4.5 Trepan																											TF_RS4.0-4.5

4.0-4.5mm Pilot Drills



4.0-4.5 Pilot 12mm																											RS_PB 4.0-4.5/12
4.0-4.5 Pilot 14mm																											RS_PB 4.0-4.5/14
4.0-4.5 Pilot 16mm																											RS_PB 4.0-4.5/16

4.0mm Cortical Drills



4.0 Cortical 3mm																											RS_KB 4.0/3
4.0 Cortical 5mm																											RS_KB 4.0/5
4.0 Cortical 7mm																											RS_KB 4.0/7

4.5mm Cortical Drills



4.5 Cortical 3mm																											RS_KB 4.5/3
4.5 Cortical 5mm																											RS_KB 4.5/5
4.5 Cortical 7mm																											RS_KB 4.5/7

4.0mm Final Drills



4.0 Final 12mm																											RS_B 4.0/12
4.0 Final 14mm																											RS_B 4.0/14
4.0 Final 16mm																											RS_B 4.0/16
4.0 Final 18mm																											RS_B 4.0/18
4.0 Final 20mm																											RS_B 4.0/20
4.0 Final 22mm																											RS_B 4.0/22
4.0 Final 24mm																											RS_B 4.0/24
4.0 Final 26mm																											RS_B 4.0/26
4.0 Final 28mm																											RS_B 4.0/28

