



1000 Corporate Drive, Marshfield, WI 54449
1-844-484-4998 · www.argondentalusa.com

Argon K3Pro™ 5.0 - 6.0mm Rapid Precision Surgical Kit

Drill Usage Chart

IMPORTANT NOTE:

To maintain optimal performance, Argon Dental USA recommends replacing drills after 20 uses. Please use this chart as a tracking device. Check the box which corresponds to the drill, each time it is used.

5.0-6.0mm Punch, Starter and Trepan



	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	Product #	
5.0-6.0 Punch																											GS_RS5.0-6.0
5.0-6.0 Starter																											RS_PF5.0-6.0
5.0-6.0 Trepan																											TF_RS5.0-6.0

5.0-6.0mm Pilot Drills



5.0-6.0 Pilot 12mm																											RS_PB 5.0-6.0/12
5.0-6.0 Pilot 14mm																											RS_PB 5.0-6.0/14
5.0-6.0 Pilot 16mm																											RS_PB 5.0-6.0/16

5.0mm Cortical Drills



5.0 Cortical 3mm																											RS_KB 5.0/3
5.0 Cortical 5mm																											RS_KB 5.0/5
5.0 Cortical 7mm																											RS_KB 5.0/7

6.0mm Cortical Drills



6.0 Cortical 3mm																											RS_KB 6.0/3
6.0 Cortical 5mm																											RS_KB 6.0/5
6.0 Cortical 7mm																											RS_KB 6.0/7

5.0mm Final Drills




5.0 Final 12mm																											RS_B 5.0/12
5.0 Final 14mm																											RS_B 5.0/14
5.0 Final 16mm																											RS_B 5.0/16
5.0 Final 18mm																											RS_B 5.0/18
5.0 Final 20mm																											RS_B 5.0/20
5.0 Final 22mm																											RS_B 5.0/22
5.0 Final 24mm																											RS_B 5.0/24

6.0mm Final Drills



	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	Product #	
6.0 Final 12mm																											RS_B 6.0/12
6.0 Final 14mm																											RS_B 6.0/14
6.0 Final 16mm																											RS_B 6.0/16
6.0 Final 18mm																											RS_B 6.0/18
6.0 Final 20mm																											RS_B 6.0/20

5.0-6.0mm Taps



5.0mm Tap																											RS_GS 5.0
6.0mm Tap																											RS_GS 6.0

5.0-6.0mm Implant Inserter
Hand Ratchets



5.0-6.0 Ratchet 3mm																											RSIE_RA 5.0-6.0/3
5.0-6.0 Ratchet 5mm																											RSIE_RA 5.0-6.0/5
5.0-6.0 Ratchet 7mm																											RSIE_RA 5.0-6.0/7

5.0-6.0mm Contra Angle
Implant Inserters



5.0-6.0 Contra 3mm																											RSIE_WS 5.0-6.0/3
5.0-6.0 Contra 5mm																											RSIE_WS 5.0-6.0/5
5.0-6.0 Contra 7mm																											RSIE_WS 5.0-6.0/7

5.0-6.0mm Implant Inserter
Hand Ratchets-Short



5.0-6.0 Ratchet-S 2.5mm																											RSIE_RA 5.0-6.0/2.5_S
5.0-6.0 Ratchet-S 4.5mm																											RSIE_RA 5.0-6.0/4.5_S
5.0-6.0 Ratchet-S 6.5mm																											RSIE_RA 5.0-6.0/6.5_S

5.0-6.0mm Contra Angle
Implant Inserters-Short



5.0-6.0 Contra-S 2.5mm																											RSIE_WS 5.0-6.0/2.5_S
5.0-6.0 Contra-S 4.5mm																											RSIE_WS 5.0-6.0/4.5_S
5.0-6.0 Contra-S 6.5mm																											RSIE_WS 5.0-6.0/6.5_S