



1000 Corporate Drive, Marshfield, WI 54449
1-844-484-4998 · www.argondentalusa.com

Argon K3Pro™ Basic Rapid Surgical Kit Drill Usage Chart

IMPORTANT NOTE:
To maintain optimal performance, Argon Dental USA recommends replacing drills after 20 uses. Please use this chart as a tracking device. Check the box which corresponds to the drill, each time it is used.

Bur, Starter, Extension & Pilot



	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	Product #	
Round Bur																											KB1.9K
Lance Drill																											KB.STLD
Drill Extension																											KBV03
2.0 Pilot Drill																											KB2.0/18.0OS

**Implant Drills
Conical Surgical Steel**



3.0 x 18.0mm Drill																											KKB3.0/18.0OS
3.5 x 18.0mm Drill																											KKB3.5/18.0OS
4.0 x 18.0mm Drill																											KKB4.0/18.0OS
4.5 x 18.0mm Drill																											KKB4.5/18.0OS
5.0 x 18.0mm Drill																											KKB5.0/18.0OS
6.0 x 18.0mm Drill																											KKB6.0/18.0OS

Implant Taps



3.0mm Tap																											GS_3.0
3.5mm Tap																											GS_3.5
4.0mm Tap																											GS_4.0
4.5mm Tap																											GS_4.5
5.0mm Tap																											GS_5.0
6.0mm Tap																											GS_6.0

Implant Sinkers



3.0mm Sinker																											K_3.0
3.5mm Sinker																											K_3.5
4.0mm Sinker																											K_4.0
4.5mm Sinker																											K_4.5
5.0mm Sinker																											K_5.0
6.0mm Sinker																											K_6.0

