

# K3Pro™ Torque Wrench Calibration Work Order Form

www.argondentalusa.com • 1-844-484-4998



## Customer Information



Date: \_\_\_\_\_

Name: \_\_\_\_\_

Name of Practice: \_\_\_\_\_

Shipping Address: \_\_\_\_\_

Phone number: \_\_\_\_\_ Email: \_\_\_\_\_

Number of Torque Wrenches Submitted with Work Order Form: \_\_\_\_\_

## Annual Torque Wrench Testing is Critical to Patient Success

Proper use and care of your torque wrench is critically important to your patient's long-term dental success. If a dental abutment is overtorqued, the screw can break inside the implant. If it is undertorqued, the abutment can loosen and work its way out.

Argon Dental USA is committed to upholding the highest standards in product manufacturing, offering annual calibration torque wrench testing.

## Calibration Testing Service

Depending on the results of the torque wrench calibration test, we will either adjust the spring tension or replace the spring to bring the torque wrench back into calibration.

We charge \$21.00 to test each torque wrench. If the results of the test show that the spring needs to be replaced, the total fee will be \$53.00 for the replacement spring and additional calibration testing. As part of our service, we will provide a report that details the initial test results and any adjustments that were made.

Please send  
this form and your  
K3Pro torque wrenches  
to the lab  
**AUTOCLAVED!**



Calibration Testing Machine



Set Screw

Spring

## Internal Use ONLY

Code for Service: ☐ PKG20 ☐ PKG50

Wrench Serial #'s: \_\_\_\_\_

Date Received: \_\_\_\_\_



Argon Dental USA™

Order To: Argon Dental USA - 1000 Corporate Drive - P.O. Box 770 - Marshfield, WI 54449